

*New, from Brasseler USA...*

# **PneuMicro™** Wiredriver

*For Precise, Dependable Control and Power*



The PneuMicro™ Wiredriver from Brasseler USA brings dependable lube-free technology to the surgeon. This pneumatically-powered unit, with a comfort grip handle design, places precise, variable speed control at your fingertips.

The PneuMicro™ Wiredriver's motor and valve delivers the control and power necessary for drilling and wire fixation, and our collet design offers expanded versatility of drilling with quick connect drills.

## *An Array of New, Outstanding Features:*

- Fatigue-reducing, ergonomic handle design.
- Lightweight and compact.
- Highly responsive, variable speed trigger.
- Improved, robust, high-performance motor design.
- Lube-free motor design for reduced maintenance.
- Expanded wire/pin diameter range...from 0.7 mm to 2.0 mm.
- Convenient, single hand directional control switch.
- Improved valve design for enhanced directional control with "positive touch."



**BRASSLER USA®**  
SURGICAL POWER & ACCESSORIES