


Standard & Synthes® Quick Connect Twist Drills

Standard Twist Drills

Manufactured from 316LVM stainless steel, our precision ground implant grade material.
Available six (6) per pack(non-sterile) and one (1) per pack(sterile).



5" Length (Non-Sterile)*

Catalog Number	Diameter - Inch (mm)	
KM166-00-00	1/16" (1.6mm)	Non-Sterile
KM166-00-01	5/64" (2.0mm)	Non-Sterile
KM166-00-02	3/32" (2.4mm)	Non-Sterile
KM166-00-03	7/64" (2.8mm)	Non-Sterile
KM166-00-04	1/8" (3.2mm)	Non-Sterile
KM166-00-05	9/64" (3.6mm)	Non-Sterile
KM166-00-06	5/32" (4.0mm)	Non-Sterile
KM166-00-07	3/16" (4.8mm)	Non-Sterile
KM166-00-08	1/4" (6.4mm)	Non-Sterile

5" Length (Sterile)

Catalog Number	Diameter - Inch (mm)	
KM166-10-00	1/16" (1.6mm)	Sterile
KM166-10-01	5/64" (2.0mm)	Sterile
KM166-10-02	3/32" (2.4mm)	Sterile
KM166-10-03	7/64" (2.8mm)	Sterile
KM166-10-04	1/8" (3.2mm)	Sterile
KM166-10-05	9/64" (3.6mm)	Sterile
KM166-10-06	5/32" (4.0mm)	Sterile
KM166-10-07	3/16" (4.8mm)	Sterile
KM166-10-08	1/4" (6.4mm)	Sterile



Order/Information - USA & Canada

1-800-535-6638 or 1-912-925-8525



BRASSLER USA
SURGICAL POWER & ACCESSORIES