

BRASSLER USA[®]

HIGH SPEED BURS



PRECISION. PERFORMANCE.

Your Value Provider





BRASSLER USA[®]
SURGICAL POWER & ACCESSORIES

Brasseler USA's high Speed burs are offered at significant cost savings and are fully compatible with the original manufacturer's equipment. Our high speed burs are designed and manufactured to very rigid specifications and are subjected to extensive engineering analysis and strenuous testing programs to ensure exceptional product performance.






Brasseler USA® High Speed Burs

Compatible with Medtronic® Midas Rex® Legend® Systems

Brasseler USA® Order Number	Replaces	Description	Head Diameter (mm)	Head Length (mm)	Brasseler USA® High Speed Burs Compatible with the following Attachments
					
BR9AC60	<i>9AC60</i>	Acorn - Fluted	6.0	7.7	AA09S, AS09, AS09S, AVA09, AVS09
BR9AC75	<i>9AC75</i>	Acorn - Fluted	7.5	9.4	AA09S, AS09, AS09S, AVA09, AVS09
BR9AC90	<i>9AC90</i>	Acorn - Fluted	9.0	11.4	AA09S, AS09, AS09S, AVA09, AVS09
BR14AC60	<i>14AC60</i>	Acorn - Fluted	6.0	7.7	AA14, AA14S, AS14, AS14S, AVA14, AVS14
BR14AC75	<i>14AC75</i>	Acorn - Fluted	7.5	9.4	AA14, AA14S, AS14, AS14S, AVA14, AVS14
					
BR7BA20	<i>7BA20</i>	Ball - Fluted	2.0	N/A	AA07, AS07 AVA07, AVS07
BR7BA50	<i>7BA50</i>	Ball - Fluted	5.0	N/A	AA07, AS07 AVA07, AVS07
BR9BA30	<i>9BA30</i>	Ball - Fluted	3.0	N/A	AA09S, AS09, AS09S, AVA09, AVS09
BR9BA40	<i>9BA40</i>	Ball - Fluted	4.0	N/A	AA09S, AS09, AS09S, AVA09, AVS09
BR9BA50	<i>9BA50</i>	Ball - Fluted	5.0	N/A	AA09S, AS09, AS09S, AVA09, AVS09
BR9BA60	<i>9BA60</i>	Ball - Fluted	6.0	N/A	AA09S, AS09, AS09S, AVA09, AVS09
BR9BA75	<i>9BA75</i>	Ball - Fluted	7.5	N/A	AA09S, AS09, AS09S, AVA09, AVS09
BR10BA30	<i>10BA30</i>	Ball - Fluted	3.0	N/A	AA10, AA10S, AS10, AS10S, AVA10, AVS10
BR10BA40	<i>10BA40</i>	Ball - Fluted	4.0	N/A	AA10, AA10S, AS10, AS10S, AVA10, AVS10
BR10BA50	<i>10BA50</i>	Ball - Fluted	5.0	N/A	AA10, AA10S, AS10, AS10S, AVA10, AVS10
BR14BA30	<i>14BA30</i>	Ball - Fluted	3.0	N/A	AA14, AA14S, AS14, AS14S, AVA14, AVS14
BR14BA40	<i>14BA40</i>	Ball - Fluted	4.0	N/A	AA14, AA14S, AS14, AS14S, AVA14, AVS14
BR14BA50	<i>14BA50</i>	Ball - Fluted	5.0	N/A	AA14, AA14S, AS14, AS14S, AVA14, AVS14




Medtronic® Midas Rex® and Legend® are registered trademarks of Medtronic, Inc.
Brasseler USA® is not affiliated with Medtronic, Inc.

Brasseler USA® High Speed Burs Compatible with Medtronic® Midas Rex® Legend® Systems

Brasseler USA® Order Number	Replaces	Description	Head Diameter (mm)	Head Length (mm)	Brasseler USA® High Speed Burs Compatible with the following Attachments
					
BR9BA50D	<i>9BA50D</i>	Ball - Diamond	5.0	N/A	AA09S, AS09, AS09S, AVA09, AVS09
BR10BA50D	<i>10BA50D</i>	Ball - Diamond	5.0	N/A	AA10, AA10S, AS10, AS10S, AVA10, AVS10
BR14BA50D	<i>14BA50D</i>	Ball - Diamond	5.0	N/A	AA14, AA14S, AS14, AS14S, AVA14, AVS14
BR9BA30DC	<i>9BA30DC</i>	Ball - Diamond Coarse	3.0	N/A	AA09S, AS09, AS09S, AVA09, AVS09
BR9BA60DC	<i>9BA60DC</i>	Ball - Diamond Coarse	6.0	N/A	AA09S, AS09, AS09S, AVA09, AVS09
BR14BA60DC	<i>14BA60DC</i>	Ball - Diamond Coarse	6.0	N/A	AA14, AA14S, AS14, AS14S, AVA14, AVS14
BR9BA60DX	<i>9BA60DX</i>	Ball - Diamond Extra Coarse	6.0	N/A	AA09S, AS09, AS09S, AVA09, AVS09
BR10BA40DX	<i>10BA40DX</i>	Ball - Diamond Extra Coarse	4.0	N/A	AA10, AA10S, AS10, AS10S, AVA10, AVS10
					
BR9CY50	<i>9CY50</i>	Cylinder - Fluted	5.0	7.9	AA09S, AS09, AS09S, AVA09, AVS09
BR10CY40	<i>10CY40</i>	Cylinder - Fluted	4.0	6.2	AA10, AA10S, AS10, AS10S, AVA10, AVS10
BR14CY50	<i>14CY50</i>	Cylinder - Fluted	5.0	7.9	AA14, AA14S, AS14, AS14S, AVA14, AVS14
BR14CY60	<i>14CY60</i>	Cylinder - Fluted	6.0	7.9	AA14, AA14S, AS14, AS14S, AVA14, AVS14
					
BR14CY65DX	<i>14CY65DX</i>	Cylinder - Diamond Extra Coarse	6.5	9.8	AA14, AA14S, AS14, AS14S, AVA14, AVS14



Brasseler USA® High Speed Burs Compatible with Medtronic® Midas Rex® Legend® Systems

Brasseler USA® Order Number	Replaces	Description	Head Diameter (mm)	Head Length (mm)	Brasseler USA® High Speed Burs Compatible with the following Attachments
					
BR8MH17	<i>8MH17</i>	Match Head - Fluted	1.7	2.4	AS08, AVS08
BR9MH30	<i>9MH30</i>	Match Head - Fluted	3.0	3.8	AA09S, AS09, AS09S, AVA09, AVS09
BR10MH22	<i>10MH22</i>	Match Head - Fluted	2.2	3.8	AA10, AA10S, AS10, AS10S, AVA10, AVS10
BR10MH30	<i>10MH30</i>	Match Head - Fluted	3.0	3.8	AA10, AA10S, AS10, AS10S, AVA10, AVS10
BR10MH30-MN	<i>10MH30-MN</i>	Match Head - Fluted (MedNEXT® Type)	3.0	3.8	AA10, AA10S, AS10, AS10S, AVA10, AVS10
BR14MH30	<i>14MH30</i>	Match Head - Fluted	3.0	3.8	AA14, AA14S, AS14, AS14S, AVA14, AVS14
BR15MH22	<i>15MH22</i>	Match Head - Fluted	2.2	3.8	AA15, AS15 AVA15, AVS15
					
BR9MH30D	<i>9MH30D</i>	Match Head - Diamond	3.0	3.8	AA09S, AS09, AS09S, AVA09, AVS09
BR10MH30D-MN	<i>10MH30D-MN</i>	Match Head - Diamond (MedNEXT® Type)	3.0	4.1	AA10, AA10S, AS10, AS10S, AVA10, AVS10
BR14MH30D	<i>14MH30D</i>	Match Head - Diamond	3.0	3.8	AA14, AA14S, AS14, AS14S, AVA14, AVS14
					
BR10OV40	<i>10OV40</i>	Oval	4.0	6.4	AA10, AA10S, AS10, AS10S, AVA10, AVS10
BR14OV40	<i>14OV40</i>	Oval	4.0	6.4	AA14, AA14S, AS14, AS14S, AVA14, AVS14



Brasseler USA® High Speed Burs Compatible with Medtronic® Midas Rex® Legend® Systems

Brasseler USA® Order Number	Replaces	Description	Head Diameter (mm)	Head Length (mm)	Brasseler USA® High Speed Burs Compatible with the following Attachments
BR8TA11	<i>8TA11</i>	Tapered	1.1	6.4	AS08, AVS08
BR14TA30	<i>14TA30</i>	Tapered	3.0	35.1	AA14, AA14S, AS14, AS14S, AVA14, AVS14
BRF1/8TA15	<i>F1/8TA15</i>	Tapered	1.5	11.1	AF01, AF01R, AS08, AVS08
BRF2/8TA23	<i>F2/8TA23</i>	Tapered	2.3	15.9	AF02, AF02R, AS08, AVS08
BR8TD114	<i>8TD114</i>	Twist Drill	1.1	4.0 Stop Depth	AS08, AVS08
BR8TD116	<i>8TD116</i>	Twist Drill	1.1	6.0 Stop Depth	AS08, AVS08
BR8TD158	<i>8TD158</i>	Twist Drill	1.5	8.0 Stop Depth	AS08, AVS08




Brasseler USA® High Speed Burs Compatible with Medtronic® Midas Rex® Legend® Systems

Brasseler USA® Order Number	Replaces	Description	Head Diameter (mm)	Head Length (mm)	Brasseler USA® High Speed Burs Compatible with the following Attachments
BRMC254	<i>MC254</i>	Metal Cutter	25.4	1.0	ASMC
BRMC30	<i>MC30</i>	Metal Cutter	3.0	18.3	ASMC

Brasseler USA's high speed burs are offered at significant cost savings and are fully compatible with the original manufacturer's equipment. Our high speed burs are designed and manufactured to very rigid specifications and are subjected to extensive engineering analysis and strenuous testing programs to ensure exceptional product performance.

Brasseler USA[®] High Speed Replacement Burs
Compatible with Stryker[®] TPS[™], Medtronic[®] Midas Rex[®]
Classic I & II and Linvatec[®] UltraPower[™] Systems







Order Number	Description	Unit	Overall Length	Head Diameter	Head Length
					
KM-M330U	Round Fluted	mm	83.0	3.0	3.0
		in	3.27	0.12	0.12
KM-M340U	Round Fluted	mm	84.0	4.0	4.0
		in	3.31	0.16	0.16
KM-M35U	Round Fluted	mm	86.0	5.1	5.1
		in	3.39	0.20	0.20
KM-M32U	Round Fluted	mm	89.0	7.6	7.6
		in	3.50	0.30	0.30
KM-M31U	Round Fluted	mm	91.0	9.1	9.1
		in	3.58	0.36	0.36
KM-AM330U	Round Fluted	mm	137.0	3.0	3.0
		in	5.39	0.12	0.12
KM-AM340U	Round Fluted	mm	137.0	4.0	4.0
		in	5.39	0.16	0.16
KM-AM35U	Round Fluted	mm	140.0	5.1	5.1
		in	5.51	0.20	0.20
KM-AM33U	Round Fluted	mm	140.0	6.0	6.0
		in	5.51	0.24	0.24
KM-AM32U	Round Fluted	mm	143.0	7.6	7.6
		in	5.63	0.30	0.30
					
KM5120-10-230	Round Fluted Extended	mm	96.0	3.0	3.0
		in	3.78	0.12	0.12
KM5120-10-240	Round Fluted Extended	mm	96.0	4.0	4.0
		in	3.78	0.16	0.16
KM5120-10-250	Round Fluted Extended	mm	93.0	5.0	5.0
		in	3.66	0.20	0.20
KM5130-10-230	Round Fluted Extended	mm	149.0	3.0	3.0
		in	5.87	0.12	0.12
KM5130-10-240	Round Fluted Extended	mm	149.0	4.0	4.0
		in	5.87	0.16	0.16
KM5130-10-250	Round Fluted Extended	mm	146.0	5.0	5.0
		in	5.75	0.20	0.20

Linvatec[®] and UltraPower[™] are trademarks of ConMed Corporation. Medtronic[®] and Midas Rex[®] are registered trademarks of Medtronic, Inc. Stryker[®] and TPS[™] are trademarks of Stryker Corporation. Brasseler USA[®] is not affiliated with any of the above.

Brasseler USA[®] High Speed Replacement Burs

Compatible with Stryker[®] TPS[™], Medtronic[®] Midas Rex[®]
Classic I & II and Linvatec[®] UltraPower[™] Systems



Order Number	Description	Unit	Overall Length	Head Diameter	Head Length
					
KM-M340DU	Round Diamond	mm	84.0	4.0	4.0
		in	3.31	0.16	0.16
					
KM-U295	Oval Aggressive Cut	mm	89.0	4.0	8.0
		in	3.50	0.16	0.31
KM5110-15-40	Oval Aggressive Cut	mm	72.0	4.0	9.5
		in	2.83	0.16	0.37
KM5130-15-40	Oval Aggressive Cut	mm	142.0	4.0	9.5
		in	5.59	0.16	0.37
					
KM-M8U	Match Head	mm	93.0	2.9	3.8
		in	3.66	0.11	0.15
KM-AM8U	Match Head	mm	147.0	2.9	3.8
		in	5.79	0.11	0.15
					
KM-M52U	Metal Cutter (carbide)	mm	95.0	3.2	17.0
	(Round End)	in	3.74	0.13	0.67

The Brasseler USA Value

Brasseler USA has set the standard for excellence in offering surgical products compatible with original manufactures power systems.

- ✓ Low Cost / High Value – Manufacturer Direct
- ✓ Single Source Provider - Product Consolidation / Standardization
- ✓ GPO / IDN Contracts
- ✓ Products Manufactured in USA
- ✓ Product Cross Referencing on website
- ✓ Manufacturer Direct National Sales Force & Inside Account Managers
- ✓ Inquire about our Free 2nd Day Express Shipping.
- ✓ Customer Service - Product Assistance Help Desk
(Extended Customer Service Hours 7:30PM EST)
- ✓ Quality Management System – ISO Certified



One Brasseler Boulevard • Savannah, Georgia 31419 • 912-925-8525
800-535-6638 • 912-921-7578 (fax) • www.BrasselerUSAMedical.com

Linvatec® and UltraPower™ are trademarks of ConMed Corporation. Medtronic®, Midas Rex® and Legend® are registered trademarks of Medtronic, Inc.
Stryker® and TPS™ are trademarks of Stryker Corporation.
Brasseler USA® is not affiliated with any of the above.